




**INSTRUMENTS FOR
LAMINECTOMY**

KERRISON




Bissbreite Width of bite Anchura de la boca Larghezza del morso	Schaftlänge Length of shaft Longitud del vástago Lunghezza dell'asta	Standardfuß Standard footplate Lámina de punzonar standard Battente standard		Dünfuß (Mikro) Thin footplate (Micro) Lámina de punzonar estrecha (Micro) Battente sottile (Micro)	
		130° (40°)	90°	130° (40°)	90°
1 mm	180 mm	240-071-01	240-051-01		
	200 mm	240-072-01	240-052-01		
2 mm	180 mm	240-071-02	240-051-02	240-071-02T	240-051-02T
	200 mm	240-072-02	240-052-02	240-072-02T	240-052-02T
	230 mm	240-073-02	240-053-02	240-073-02T	240-053-02T
3 mm	180 mm	240-071-03	240-051-03	240-071-03T	240-051-03T
	200 mm	240-072-03	240-052-03	240-072-03T	240-052-03T
	230 mm	240-073-03	240-053-03	240-073-03T	240-053-03T
4 mm	180 mm	240-071-04	240-051-04	240-071-04T	240-051-04T
	200 mm	240-072-04	240-052-04	240-072-04T	240-052-04T
	230 mm	240-073-04	240-053-04	240-073-04T	240-053-04T
5 mm	180 mm	240-071-05	240-051-05	240-071-05T	240-051-05T
	200 mm	240-072-05	240-052-05	240-072-05T	240-052-05T
	230 mm	240-073-05	240-053-05	240-073-05T	240-053-05T
6 mm	180 mm	240-071-06	240-051-06	240-071-06T	240-051-06T
	200 mm	240-072-06	240-052-06	240-072-06T	240-052-06T





KERRISON




Bissbreite Width of bite Anchura de la boca Larghezza del morso	Schaftlänge Length of shaft Longitud del vástago Lunghhezza dell'asta	Standardfuß Standard footplate Lámina de punzonar standard Battente standard	
		130° (40°)	90°
2 mm	180 mm	240-081-02	240-061-02
	200 mm		240-062-02
3 mm	180 mm	240-081-03	240-061-03
	200 mm		240-062-03
4 mm	180 mm	240-081-04	240-061-04
	200 mm		240-062-04
5 mm	180 mm	240-081-05	240-061-05
	200 mm		240-062-05
6 mm	180 mm		240-061-06
	200 mm		240-062-06

KERRISON




Bissbreite Width of bite Anchura de la boca Larghezza del morso	Schaftlänge Length of shaft Longitud del vástago Lunghezza dell'asta	Standardfuß Standard footplate Lámina de punzonar standard Battente standard		Dünnfuß (Mikro) Thin footplate (Micro) Lámina de punzonar estrecha (Micro) Battente sottile (Micro)	
		130° (40°)	90°	130° (40°)	90°
					
1 mm	180 mm	249-071-01	249-051-01		
	200 mm	249-072-01	249-052-01		
2 mm	180 mm	249-071-02	249-051-02	249-071-02T	249-051-02T
	200 mm	249-072-02	249-052-02	249-072-02T	249-052-02T
3 mm	180 mm	249-071-03	249-051-03	249-071-03T	249-051-03T
	200 mm	249-072-03	249-052-03	249-072-03T	249-052-03T
4 mm	180 mm	249-071-04	249-051-04	249-071-04T	249-051-04T
	200 mm	249-072-04	249-052-04	249-072-04T	249-052-04T
5 mm	180 mm	249-071-05	249-051-05	249-071-05T	249-051-05T
	200 mm	249-072-05	249-052-05	249-072-05T	249-052-05T
6 mm	180 mm	249-071-06	249-051-06	249-071-06T	249-051-06T
	200 mm	249-072-06	249-052-06	249-072-06T	249-052-06T

KERRISON




Bissbreite Width of bite Anchura de la boca Larghezza del morso	Schaftlänge Length of shaft Longitud del vástago Lunghezza dell'asta	Standardfuß Standard footplate Lámina de punzonar standard Battente standard	Dünnfuß (Mikro) Thin footplate (Micro) Lámina de punzonar estrecha (Micro) Battente sottile (Micro)	mit Keramikbeschichtung / with ceramic coating con una capa de cerámica / con una copertura in cerámica	
		Standardfuß Standard footplate Lámina de punzonar standard Battente standard	Dünnfuß (Mikro) Thin footplate (Micro) Lámina de punzonar estrecha (Micro) Battente sottile (Micro)	Standardfuß Standard footplate Lámina de punzonar standard Battente standard	Dünnfuß (Mikro) Thin footplate (Micro) Lámina de punzonar estrecha (Micro) Battente sottile (Micro)
		130° (40°)	130° (40°)	130° (40°)	130° (40°)
2 mm	180 mm	240-271-02	240-271-02T	240-271-02CC	240-271-02TCC
	200 mm	240-272-02	240-272-02T	240-272-02CC	240-272-02TCC
3 mm	180 mm	240-271-03	240-271-03T	240-271-03CC	240-271-03TCC
	200 mm	240-272-03	240-272-03T	240-272-03CC	240-272-03TCC
4 mm	180 mm	240-271-04	240-271-04T	240-271-04CC	240-271-04TCC
	200 mm	240-272-04	240-272-04T	240-272-04CC	240-272-04TCC
5 mm	180 mm	240-271-05	240-271-05T	240-271-05CC	240-271-05TCC
	200 mm	240-272-05	240-272-05T	240-272-05CC	240-272-05TCC








LOVE-KERRISON




Bissbreite Width of bite Anchura de la boca Larghezza del morso	Schaftlänge Length of shaft Longitud del vástago Lunghezza dell'asta	mit Keramikbeschichtung / with ceramic coating con una capa de cerámica / con una copertura in cerámica			
		Standardfuß Standard footplate Lámina de punzonar standard Battente standard	Dünnfuß (Mikro) Thin footplate (Micro) Lámina de punzonar estrecha (Micro) Battente sottile (Micro)	Standardfuß Standard footplate Lámina de punzonar standard Battente standard	Dünnfuß (Mikro) Thin footplate (Micro) Lámina de punzonar estrecha (Micro) Battente sottile (Micro)
		130° (40°)	130° (40°)	130° (40°)	130° (40°)
2 mm	160 mm	240-090-02	240-090-02T	240-090-02CC	240-090-02TCC
3 mm	160 mm	240-090-03	240-090-03T	240-090-03CC	240-090-03TCC
4 mm	160 mm	240-090-04	240-090-04T	240-090-04CC	240-090-04TCC
5 mm	160 mm	240-090-05	240-090-05T	240-090-05CC	240-090-05TCC



mit Keramikbeschichtung / with ceramic coating
 con una capa de cerámica / con una cobertura in ceràmica

					
	2,0 x 8,0 mm	240-590-02	240-591-02	240-590-02CC	240-591-02CC
	3,0 x 9,0 mm	240-590-03	240-591-03	240-590-03CC	240-591-03CC
	4,0 x 10,0 mm	240-590-04	240-591-04	240-590-04CC	240-591-04CC


CUSHING



240-550-02



240-600-02



240-650-02

a = 130 mm

a = 150 mm

a = 180 mm




240-551-02



240-601-02



240-651-02



240-552-02




240-602-02



240-652-02


LOVE-GRUENWALD



240-550-03



240-600-03




240-650-03


a = 130 mm

a = 150 mm

a = 180 mm



240-551-03



240-601-03



240-651-03



240-552-03




240-602-03




240-652-03

SPURLING








a = 130 mm	240-550-04	240-600-04	240-650-04
a = 150 mm	240-551-04	240-601-04	240-651-04
a = 180 mm	240-552-04	240-602-04	240-652-04

CASPAR



Maulteile gefenestert und gezahnt
Jaws fenestrated and serrated
Bocas fenestradas y dentadas
Morsi fenestrati e dentellati



 2,0 x 12,0 mm	a = 140 mm	240-509-02		
	a = 160 mm	240-510-02		
	a = 185 mm	240-511-02		
 3,0 x 12,0 mm	a = 140 mm	240-509-03	240-513-03	240-516-03
	a = 160 mm	240-510-03	240-514-03	240-517-03
	a = 185 mm	240-511-03	240-515-03	240-518-03
 4,0 x 14,0 mm	a = 140 mm	240-509-04	240-513-04	240-516-04
	a = 160 mm	240-510-04	240-514-04	240-517-04
	a = 185 mm	240-511-04	240-515-04	240-518-04
 5,0 x 14,0 mm	a = 160 mm	240-510-05	240-514-05	
	a = 185 mm	240-511-05	240-515-05	
 6,0 x 16,0 mm	a = 160 mm	240-510-06		
	a = 185 mm	240-511-06		

TONTARRA
MEDIZINTECHNIK GMBH

Kantstraße 33
D-78573 Wurmlingen
Germany

Telefon: ++ 74 61 / 9 65 76 - 0
Telefax: ++ 74 61 / 9 65 76 26

e-mail: info@tontarra.de
internet: www.tontarra.de