


## The "cleanwave" System:

Especially the wavelike shape of the separating plane allows finding impurities and stains on the sliding shaft instrument which can be removed and cleaned when being detected.



The effectiveness of a standardized cleaning procedure of the cleanwave instruments was evaluated and documented by an independent hygiene laboratory.

The excellent cleaning results can be attributed to the fact that any type of cleaning utensil reaches very easily the area between upper and lower part of the sliding shaft instrument.

Beyond the cleaning aspect thanks to the wavelike shape there is a clear reduction of friction of the sliding parts which allows the surgeon a very precise and consistent procedure.

clean  wave