

GO JACK

Det første cirkulære Medical Device
i en række af bæredygtige innovationer



Holder instrumenter åbne under vaskeprocessen

Designet som en "jack" for at holde hængslede instrumenter åbne under vaskeprocessen.

Fremstillet af recirkuleret pakkemateriale

Artikel nummer: GC-01-01
Pakningsstørrelse: 30 stk.



En pilotundersøgelse har vist, at det er muligt at genanvende og oparbejde affald, såsom brugt pakkemateriale, en termoplastisk polymer, til et råmateriale af høj kvalitet til sprøjtestøbning af medicinsk udstyr (data på fil).

opitek

OUR STORY - Time to make a difference!


Healthcare is responsible for a large amount of waste and CO₂ emissions. Approximately 30% of all hospital waste are plastics and approximately 20% comes from the Operating Room. The hospitals in The Netherlands alone contribute 7% of the total national CO₂ emissions.

A disposable society, what a waste!

In 2019 Van Straten Medical co-founded the company GreenCycl with the mission to support healthcare, specifically hospitals and the medical industry, to become more sustainable.

Through a circular process, GreenCycl collects & reprocesses hospital waste via urban mining. Van Straten Medical manufactures and supplies the recycled & circular medical products.

CIRCULAR PROCESS



For our circular products, we are using sustainable packaging.
A vast amount of circular products are currently being developed.

Manufactured by Van Straten Medical. Endorsed by GreenCycl

VAN STRATEN MEDICAL

'Providing Value to Life'